Kinesiologie-Tape richtig anwenden

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support@hula-hoop-shop.eu
<https://www.hula-hoop-shop.eu>

Using Kinesiology Tape Correctly

Preparation, safety and 20 illustrated taping examples

Hoopmania Customer Instructions

These customer instructions are intended for people who have purchased kinesiology tape and want to understand basic safety rules and example tape paths. Kinesiology tape is also commonly called kinesiotape.

- Prepare clean, dry and unbroken skin without cream or oil; cut the tape and round the corners.
- Test a small area first if your skin is sensitive. Apply anchors and ends without stretch; do not stretch the tape to its limit or wrap it tightly.
- Rub the tape gently. To remove it, moisten it and peel it off slowly in the direction of hair growth while supporting the skin.
- Remove it immediately if you notice burning, itching, increasing redness, pain, numbness, tingling, coldness, discolouration or swelling.

Important information

These instructions show possible tape cuts and paths. They do not replace diagnosis, treatment or personal instruction. Seek medical or physiotherapy advice first for acute injuries, unexplained or persistent symptoms, sudden swelling or after a procedure.

Quick guide: neck, shoulder and back (1-6) · arm and hand (7-14) · abdomen, leg, knee and ankle (15-20)

1. Neck area

Guidance: basic orientation - have position and stretch demonstrated first

Required: a Y-shaped strip and a slightly shorter I-strip

The image shows the original path: tilt the head forwards, place the base of the Y-strip and guide both tails along the neck area. The shorter cross strip lies across the upper back. Apply the ends without stretch; have the position and amount of stretch demonstrated before first use.

2. Shoulder area


Guidance: basic orientation - use with care

Required: a Y-shaped tape strip

The image shows the base below the shoulder area and the two tails around the shoulder. Move the arm backwards and then forwards during application. Apply the ends without stretch. Agree the exact position and stretch with a qualified professional before first use.

3. Back area


Guidance: partner application or after instruction

Required: one I-strip measured from the shoulder area to the iliac crest

Lean the upper body slightly forwards and raise the arm over the head. The images show the strip running from the lower back along the side of the back. Seek professional advice first if pain, marked restriction of movement or an unexplained symptom is present.

4. Front chest area

Guidance: only after practical instruction

Required: a Y-shaped tape strip

The images show the original position across the front chest and shoulder area. The arm is raised to the side and moved outwards for application. Because this is a sensitive body area, have the exact technique demonstrated and do not derive it from the photographs alone.

5. Thoracic spine area

Guidance: partner application or after instruction

Required: three I-strips

The image shows two vertical strips with a cross strip over them. The earlier general instruction to use strong or maximum stretch is not adopted. Body position, path and stretch must suit the selected technique and the individual; apply the tape ends without stretch.

6. Lower back area

Guidance: partner application or after instruction

Required: four I-strips

The image shows four crossed strips on the lower back. The upper body is bent forwards for this partner application. Because of skin tension, overlaps and the hard-to-see position, have the technique demonstrated. Do not use a general maximum-stretch rule or self-treat acute back pain.

7. Fan application on the arm

Guidance: only after professional instruction

Required: two fan-cut tape strips

The image shows one fan on the upper arm and another below the elbow. Fan applications are often shown in connection with swelling. Do not tape unexplained or sudden swelling yourself; have the cause assessed medically and let qualified personnel determine path, direction and stretch.

8. Wrist area

Guidance: basic orientation - do not constrict

Required: one short I-strip

The images show a strip across the palm and wrist area. Tingling, numbness or a hand that falls asleep can have different causes. Do not use maximum stretch and seek medical or professional advice before applying tape for these symptoms.

9. Front of the upper arm

Guidance: basic orientation - agree the path with a professional

Required: two tape tails along the front of the upper arm

Apply the tape tails along the arm as shown. Place the ends without stretch. Arm position, tape path and stretch depend on the professionally selected application technique.

10. Outer forearm and elbow


Guidance: use with care - agree the stretch with a professional

Required: two I-strips

The images show a long strip from the back of the hand along the outer forearm and a short cross strip below the elbow. Bend the wrist down for application. Do not infer the stretch from the image alone; have persistent elbow symptoms assessed first.

11. Inner forearm and elbow

Guidance: use with care - monitor circulation

Required: one longer I-strip and one cross strip about 15 cm long

Position the strips as shown. The short strip lies at the wrist but must not encircle it tightly or constrict it. Remove the tape immediately if pain, tingling, numbness, coldness or discolouration occurs.

12. Thumb area

Guidance: advanced - have it demonstrated first

Required: one Y-shaped strip and one I-strip about 15 cm long

The three-part image sequence shows the cut and path of a more complex application around the thumb and wrist. Have the technique demonstrated step by step before first use. The tape must not restrict movement or circulation.

13. Back of the upper arm

Guidance: advanced - partner application

Required: a Y-shaped strip with an opening for the elbow

The images show the base at the wrist, the opening over the elbow and both tails along the back of the upper arm. Have arm position, cut and path demonstrated in practice. The opening and tape tails must not constrict; apply the ends without stretch.

14. Finger area

Guidance: only after professional instruction

Required: three narrow tape strips, about half the standard tape width

The images show a lengthwise strip from the wrist to the finger and two short securing strips. If a capsule, ligament or bone injury is suspected, do not self-treat from the photograph; seek medical assessment. The application must not constrict or impair circulation.

15. Abdominal area - path example only, with no efficacy claim


Guidance: path example only - no efficacy claim

Required: one I-strip measured between the iliac crest and rib margin

This example shows only one possible tape path on the abdominal surface. It does not imply an effect on digestion, internal organs, pain or other symptoms. Abdominal symptoms, constipation or other persistent complaints require medical assessment.

16. Fan application on the leg - only after professional instruction

Guidance: only after professional instruction

Required: two fan-cut strips measured on the leg

The images show one fan from the upper leg to the knee and another from the knee towards the heel. Swelling or fluid retention must first be assessed to determine the cause. Qualified personnel should determine the path, direction and stretch of this application.

17. Knee area

Guidance: only after professional instruction

Required: four I-strips

The image sequence shows several strips around the bent knee. The earlier general recommendation to use high stretch is not adopted. Knee position, order, path and stretch must be determined professionally. Do not self-tape an acute injury, instability, marked swelling or persistent pain.

18. Kneecap area - example path after instruction

Guidance: only after diagnosis and professional instruction

Required: a Y-shaped tape strip

The image shows the base on the outside of the knee and two tails running beside the kneecap. This is not an instruction for treating an injury. Use this type of application only after diagnosis and professional instruction.

19. Thigh area

Guidance: only after professional instruction

Required: four short I-strips, each about 15 cm long

The image shows four strips in a crossed arrangement. The earlier general stretch instruction is not adopted. If a strain or another acute injury is suspected, first have the cause and safe loading assessed; stretch and position follow professional guidance and skin tolerance.

20. Ankle area - multi-part application after instruction

Guidance: only after professional instruction

Required: four I-strips, three of them shorter

The images show a multi-part ankle application. They do not replace an examination after a twist, fall or instability. Cut, foot position, order, path and stretch must be demonstrated professionally. Remove the tape immediately if pain, numbness, tingling, coldness or discolouration increases.

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Hoopomania GmbH