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Hoopomania Rückenstrecker

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ENGLISH USER MANUAL

Hoopomania Back Stretcher

Assembly, adjustment and controlled use

RS-300 · with magnetic elements


RS-200 · with massage nodules

RS-400 · with acupressure nodules

Important: The back stretcher is not a substitute for medical assessment or treatment. Do not use it in the event of acute or unexplained pain, injury, numbness, signs of paralysis, dizziness, neurological symptoms, after surgery or when uncertain; seek professional advice.

1. Product and height settings

Base and support arch

The back stretcher consists of a base and a flexible support arch. The three slots in the base change the arch height. Always begin at level 1 and adjust the device only when it is not bearing any weight.

The surface and contact points may vary depending on the model supplied.

Transport position
approx. 40 × 25 × 3 cm

Level 1
Starting position, approx. 7 cm

Level 2
Higher position, approx. 9 cm

Level 3
Highest position, approx. 11 cm

Support arch

The support arch carries the backrest. Check the top and underside on the model supplied and do not fit it the wrong way round.

Orientation

The narrow end points towards the head. The wide end sits at the lower back or pelvic end.

Levels 2 and 3: Use them only when assembly, lying down and getting up at level 1 can be completed calmly and without pain or significant discomfort. A higher level is not required.

2. Assembly

Place the base and support arch on a stable, level, non-slip surface with the wide end towards you. The back stretcher must not bear any weight during assembly.

1. Align the parts

Pre-bend the support arch slightly. Align the wide ends of the base and support arch with each other.

2. Insert the wide end

Hold the base steady with your knees. Slide the wide end of the support arch into the single slot.

3. Position your hands

Place your thumbs in the opening; your fingers rest on top of the support arch. Guide the part evenly and do not twist it.

4. Lock into level 1

Press the narrow end down in a controlled manner and slide it gently into the level 1 slots. Both sides must be seated fully and evenly.

Check after assembly

- The support arch is secure at both ends and set to the same level.
- The device sits firmly on the surface without rocking; there are no cracks, sharp edges, loose parts or significant deformation.
- The narrow end will point towards the head and the wide end towards the lower back.
- Change the level only when the device is not bearing any weight; never adjust it while lying on it.

Do not use if the support arch does not lock securely, parts are damaged or deformed, or the device slips on the surface.

Assembly complete

Check the selected level and orientation once more. Only then continue with the guidance.

3. Guidance for use

Read and follow this guidance in full before using the back stretcher.

Before use

- Use only for the controlled positioning described in this manual.
- Use on a level, non-slip, load-bearing surface; do not use on a bed or soft sofa.
- Familiarise yourself with all illustrations, positions and the safe way to sit up before use.
- Seek professional or medical advice before use if you have an existing back or neck injury, a known spinal condition, osteoporosis, severe or unexplained pain, neurological symptoms, have recently had surgery, or are unsure whether use is appropriate for you.
- Begin at level 1. Select level 2 or 3 only if the lower position is secure, calm and possible without pain or significant discomfort.
- Breathe calmly, allow your body to rest on the device without additional pressure, and do not force any position or stretch.
- The narrow end points towards the head; the wide end points towards the lower back or pelvis.
- Stop immediately if you experience pain, tingling, numbness, dizziness, breathing problems or marked discomfort.

Duration and frequency

There is no fixed duration suitable for everyone. Start briefly and remain in a position only while it is calm and possible without pain or significant discomfort. Longer or more frequent use is not automatically better. Stop whenever any warning sign occurs.

4. Using the back stretcher

These positions are options, not a compulsory sequence. Use only a controlled range that does not cause discomfort.


1. Sit down and lower yourself slowly

Sit in front of the arch with your knees bent and support yourself at the sides. Lower your upper body slowly and centre your back; the wide end sits at the lower back. Place a thin pillow under your head if needed.

2. Basic position

Keep your knees bent, rest your arms at your sides and breathe calmly. This low position is the starting point. Continue only if it remains stable and free from pain or significant discomfort.

3. Moderate position

Move your arms over your head and shoulders only as far as you can without pain or significant discomfort. Your elbows do not need to reach the floor. Avoid jerky movements.

4. Extended position

Extend your arms and legs further only if the previous position is secure and comfortable. A pillow may be removed, but this is not required. Do not force a maximum stretch.

5. Roll onto your side and sit up

To finish, roll onto your side, remove the device and sit up in stages. Remain seated briefly and stand only when you are not dizzy or unwell.

Stopping, cleaning and storage

- Stop using the device and seek professional advice if unusual or persistent symptoms occur.
- Clean with a soft, slightly damp cloth and allow to dry completely.
- Store dry, without any load and out of reach of children; do not leave permanently under tension.
- Do not continue using a device that is damaged, cracked, significantly deformed or does not lock securely.

Final note

The back stretcher is intended for controlled positioning and movement as described in this manual. No claim of prevention, healing or pain relief is made. Use only a setting that is possible without pain or significant discomfort, and put the device away after use.

Manufacturer
Hoopomania GmbH
Käthe-Paulus-Str. 1
85092 Kösching · Germany

Contact
Email: support@hula-hoop-shop.eu
Web: <https://www.hula-hoop-shop.eu>

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